

For Research Use Only

EPHA1 Monoclonal antibody

Catalog Number: 65024-1-Ig



Basic Information

Catalog Number:

65024-1-Ig

Size:

150ul, Concentration: 5000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC130291

GeneID (NCBI):

2041

UNIPROT ID:

P21709

Full Name:

EPH receptor A1

Calculated MW:

976 aa, 108 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:20-1:200

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : Hela cells,

IHC : Breast cancer tissue,

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

