

For Research Use Only

Anti-Mouse CD115 (c-fms) (AFS98)

Catalog Number: 65095-1-Ig



Basic Information

Catalog Number: 65095-1-Ig	GenBank Accession Number: BC036343	Purification Method: Affinity purification
Size: 100ug, 500 µg/ml	GeneID (NCBI): 12978	CloneNo.: AFS98
Source: Rat	Full Name: colony stimulating factor 1 receptor	
Isotype: IgG2a, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide, pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug Anti-Mouse CD115 (c-fms) (65095-1-Ig, Clone: AFS98) and FITC anti-rat IgG2a Antibody at dilution 1:200 (green), or stained with 0.5 ug isotype control and FITC anti-rat IgG2a Antibody at dilution 1:200 (black). Cells were not fixed.