

For Research Use Only

Anti-Mouse LPAM-1 (DATK32)

Catalog Number: 65101-1-Ig



Basic Information

Catalog Number: 65101-1-Ig	GenBank Accession Number: BC011184	Purification Method: Affinity purification
Size: 100ug, 7.6 mg/ml	GeneID (NCBI): 16421	CloneNo.: DATK32
Source: Rat	Full Name: integrin beta 7	
Isotype: Rat IgG2a, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at 2-8°C.

Storage Buffer:
Unkonwn

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data