

For Research Use Only

Anti-Mouse LPAM-1 (DATK32)

Catalog Number: 65101-1-Ig



Basic Information

Catalog Number:

65101-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC011184

GeneID (NCBI):

16421

Full Name:

integrin beta 7

Purification Method:

Affinity purification

CloneNo.:

DATK32

Applications

Tested Applications:

FC

Species Specificity:

mouse

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse bone marrow cells were surface stained with 1 ug Anti-Mouse LPAM-1 (DATK32) (65101-1-Ig, Clone: DATK32) or 1 ug Rat IgG2a Isotype Control (2A3) (65209-1-Ig, Clone: 2A3) and FITC anti-Rat IgG2a Antibody. Cells were co-stained with APC Anti-Mouse CD45R (B220). Cells were not fixed.