

For Research Use Only

Anti-Mouse CD275 (B7-H2) (HK5.3)

Catalog Number: 65114-1-Ig



Basic Information

Catalog Number: 65114-1-Ig	GenBank Accession Number: BC029227	Purification Method: Affinity purification
Size: 100ug,	GeneID (NCBI): 50723	CloneNo.: HK5.3
Source: Rat	Full Name: icos ligand	
Isotype: Rat IgG2a, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at 2-8°C.
Storage Buffer:
Unkonwn

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data