

For Research Use Only

# Anti-Mouse CD24 (M1/69)

Catalog Number: 65121-1-Ig



## Basic Information

Catalog Number:

65121-1-Ig

Size:

100ug , 500 µg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC075622

GeneID (NCBI):

12484

Full Name:

CD24a antigen

Purification Method:

Affinity purification

CloneNo.:

M1/69

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (M1/69) (65121-1-Ig, Clone: M1/69) or 0.5 ug Rat IgG2b Isotype Control (LTF-2) (65211-1-Ig, Clone: LTF-2), and FITC anti-rat IgG2b Antibody. Cells were not fixed.