

For Research Use Only

Anti-Human IgD (IA6-2)

Catalog Number: 65175-1-Ig



Basic Information

Catalog Number:

65175-1-Ig

Size:

100ug, 1 mg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC021276

GeneID (NCBI):

3495

Full Name:

immunoglobulin heavy constant

delta

Calculated MW:

573 aa, 63 kDa

Purification Method:

Purified by protein-A affinity chromatography

CloneNo.:

IA6-2

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), 1 ug Anti-Human IgD (65175-1-Ig, Clone: IA6-2) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:500. Cells were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), 1 ug Mouse IgG2a isotype control and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:500. Cells were not fixed.