

For Research Use Only

# Anti-Mouse CD41 (MWReg30)

Catalog Number: 65189-1-Ig



## Basic Information

Catalog Number:

65189-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120493

GeneID (NCBI):

16399

UNIPROT ID:

Q9QUM0

Full Name:

integrin alpha 2b

Purification Method:

Purified by protein-G affinity chromatography

CloneNo:

MWReg30

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> C57BL/6 mouse peripheral blood platelets were surface stained with 0.5 ug Anti-Mouse CD41 (65189-1-Ig, Clone: MWReg30) and FITC anti-rat IgG1 Antibody (red), or 0.5 ug Rat IgG1 Isotype Control (HRPN) (65212-1-Ig, Clone: HRPN) and FITC anti-rat IgG1 Antibody (blue). Cells were not fixed.