

For Research Use Only

# Anti-Mouse CD48 (HM48-1)

Catalog Number: 65227-1-Ig



---

## Basic Information

<b>Catalog Number:</b> 65227-1-Ig	<b>GenBank Accession Number:</b> BC060977	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 500ug, 0.5 mg/mL	<b>GeneID (NCBI):</b> 12506	<b>CloneNo.:</b> HM48-1
<b>Source:</b> Armenian Hamster	<b>Full Name:</b> CD48 antigen	
<b>Isotype:</b> IgG		

---

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Mouse

---

## Background Information

---

### Storage

**Storage:**  
Store at 2-8°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.

---

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD48 (65227-1-Ig, Clone:HM48-1) and FITC anti-Armenian Hamster IgG Antibody at dilution 1:100 (red), or 0.5 ug Control Antibody. Cells were not fixed.