

For Research Use Only

# Anti-Rat CD43 (W3/13)

Catalog Number: 65237-1-Ig



## Basic Information

Catalog Number:

65237-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

24796

UNIPROT ID:

P13838

Full Name:

sialophorin

Calculated MW:

38 kDa

Purification Method:

Affinity purification

CloneNo.:

W3/13

## Applications

Tested Applications:

FC

Species Specificity:

Rat

## Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> Wistar rat splenocytes were surface stained with 1 ug Anti-Rat CD43 (65237-1-Ig, Clone:W3/13) and FITC anti-Mouse IgG1 Antibody at dilution 1:100, or 1 ug Control Antibody. Cells were not fixed.