

For Research Use Only

Anti-Mouse CD94 (18D3)

Catalog Number: 65240-1-Ig



Basic Information

Catalog Number:

65240-1-Ig

Size:

100ug, 0.5 mg/mL

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,
subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface co-stained with APC Anti-Mouse NK1.1 (CD161) (PK136) and 0.5 ug Anti-Mouse CD94 (65240-1-Ig, Clone:18D3) or 0.5 ug Rat IgG2a Isotype Control (2A3), and FITC anti-Rat IgG2a Antibody at dilution 1:100. Cells were not fixed.