

For Research Use Only

# Anti-Human CD155 (SKII.4)

Catalog Number: 65251-1-Ig



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
65251-1-Ig	BC015542	Affinity purification
Size:	GeneID (NCBI):	CloneNo.:
100ug, 0.5 mg/mL	5817	SKII.4
Source:	Full Name:	
Mouse	poliovirus receptor	
Isotype:	Calculated MW:	
IgG1, kappa	45 kDa	

## Applications

Tested Applications:  
FC, ELISA

Species Specificity:  
Human

## Background Information

### Storage


Storage:  
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> U937 cells were surface stained with 0.2 ug Anti-Human CD155 (65251-1-Ig, Clone:SKII.4) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed.