

For Research Use Only

Anti-Mouse GL7 (GL-7)

Catalog Number: 65261-1-Ig



Basic Information

Catalog Number: 65261-1-Ig	GenBank Accession Number: GeneID (NCBI): 104230	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	Full Name: lymphocyte antigen 77	CloneNo.: GL-7
Source: Rat		
Isotype: IgM, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Human, Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ unstimulated (blue) or plate-bound anti-CD3 and anti-CD28 stimulated (red) mouse splenocytes were surface stained with 0.5 μ g Anti-Mouse GL7 (65261-1-Ig, Clone:GL-7) and FITC anti-rat IgM Antibody at dilution 1:100. Cells were not fixed.