

For Research Use Only

# Anti-Human CD55 (JS11)

Catalog Number: 65271-1-Ig



## Basic Information

<b>Catalog Number:</b> 65271-1-Ig	<b>GenBank Accession Number:</b> BC001288	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100ug, 0.5 mg/ml	<b>GeneID (NCBI):</b> 1604	<b>CloneNo.:</b> JS11
<b>Source:</b> Mouse	<b>Full Name:</b> CD55 molecule, decay accelerating factor for complement (Cromer blood group)	
<b>Isotype:</b> IgG1	<b>Calculated MW:</b> 41 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

### Storage


**Storage:**  
Store at 2-8°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug Anti-Human CD55 (65271-1-Ig, Clone:JS11K) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed. Lymphocytes were gated.