

For Research Use Only

Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: 65571-1-MR



Basic Information

Catalog Number:

65571-1-MR

Size:

100ug, 1 mg/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Applications

Tested Applications:

FC

Species Specificity:

n/a

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with 0.2 μ g Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody (65571-1-MR, Clone:MOPC-173) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.2 μ g Mouse IgG2a Isotype Control (MOPC-173) (65244-1-Ig, Clone: MOPC-173) (blue), or 0.2 μ g Anti-Human CD8 (UCHT4) (65204-1-Ig, Clone: UCHT4)(green). Cells were not fixed. Lymphocytes were gated.