

For Research Use Only

Anti-Human STRO-1 (STRO-1) Mouse IgG2a Recombinant Antibody

Catalog Number: 65585-1-MR



Basic Information

Catalog Number:

65585-1-MR

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Protein A purification

Size:

100ug, 1 mg/ml

Full Name:

CloneNo.:

STRO-1

Source:

Mouse

Isotype:

IgG2a

Applications

Tested Applications:

FC

Species Specificity:

human

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human blood cells were surface unstained or stained with 0.5 ug Anti-Human STRO-1 Mouse IgG2a Recombinant Antibody (65585-1-MR, Clone: STRO-1) and CoraLite®488-Conjugated AffiniPure Donkey Anti-Mouse IgG(H+L). Cells were not fixed. CD45 dim cells were gated.



1x10⁶ human blood cells were unstained (blue) or surface stained with 0.5 ug Anti-Human STRO-1 Mouse IgG2a Recombinant Antibody (65585-1-MR, Clone: STRO-1) and CoraLite®488-Conjugated AffiniPure Donkey Anti-Mouse IgG(H+L) (red). Cells were not fixed. CD45 dim cells were gated.