

For Research Use Only

# Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody

Catalog Number: 65626-1-MR



## Basic Information

Catalog Number:

65626-1-MR

Size:

100ug, 1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

Full Name:

CD1c molecule

Purification Method:

Protein A purification

CloneNo.:

L161

## Applications

Tested Applications:

FC

Species Specificity:

human

## Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human peripheral blood lymphocytes were surface stained with 0.25 ug Anti-Human CD1c (L161) Mouse IgG2a Recombinant Antibody (65626-1-MR, Clone: L161) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) and PE-conjugated Goat Anti-Mouse IgG. Cells were co-stained with CoraLite® Plus 647 Anti-Human CD19. Cells were not fixed.