

For Research Use Only

# Anti-Mouse PD-L1/CD274 (10F.9G2) Mouse IgG2a Recombinant Antibody, PBS Only

Catalog Number: 65685-1-PBS



## Basic Information

Catalog Number:

65685-1-PBS

Size:

1mg, 3 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC066841

GeneID (NCBI):

60533

Full Name:

CD274 antigen

Purification Method:

Protein A purification

CloneNo.:

10F.9G2

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10<sup>6</sup> untreated or anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-L1/CD274 (10F.9G2) Mouse IgG2a RecAb (65685-1-MR, Clone:10F.9G2) and PE-conjugated Goat Anti-Mouse IgG. Cells were co-stained with APC Anti-Mouse CD3. Cells were not fixed. This data was developed using the same antibody clone with 65685-1-PBS in a different storage buffer formulation.

1x10<sup>6</sup> anti-CD3/CD28 treated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-L1/CD274 (10F.9G2) Mouse IgG2a RecAb (65685-1-MR, Clone: 10F.9G2) (red) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) (blue) and PE-conjugated Goat Anti-Mouse IgG. 1x10<sup>6</sup> untreated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse PD-L1/CD274 (10F.9G2) Mouse IgG2a RecAb (65685-1-MR,

