

For Research Use Only

TST Monoclonal antibody

Catalog Number: 66018-2-Ig



Basic Information

Catalog Number: 66018-2-Ig	GenBank Accession Number: BC010148	Purification Method: Protein A purification
Size: 150ul , Concentration: 447 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 7263	CloneNo.: 1D3D7
Source: Mouse	UNIPROT ID: Q16762	
Isotype: IgG2b	Full Name: thiosulfate sulfurtransferase (rhodanese)	
Immunogen Catalog Number: AG9235	Calculated MW: 297 aa, 33 kDa	
	Observed MW: 30-33 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data