

For Research Use Only

PI4KA Monoclonal antibody

Catalog Number: 66027-1-Ig



Basic Information

Catalog Number: 66027-1-Ig	GenBank Accession Number: BC018120	Purification Method: Protein A purification
Size: 150ul	GeneID (NCBI): 5297	CloneNo.: 4E8C11
Source: Mouse	Full Name: phosphatidylinositol 4-kinase, catalytic, alpha	
Isotype: IgG2b	Calculated MW: 2044 aa, 231 kDa	
Immunogen Catalog Number: AG17168		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

PI4KA, also named as PIK4, PIK4CA, belongs to the PI3/PI4-kinase family and Type III PI4K subfamily. It acts on phosphatidylinositol (PtdIns) in the first committed step in the production of the second messenger inositol-1,4,5-trisphosphate. PI4KA plays a role in the formation of membrane complexes where HCV replication takes place. (PMID:19605471).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data