

For Research Use Only

PI4KA Monoclonal antibody

Catalog Number: 66027-1-Ig



Basic Information

Catalog Number:	66027-1-Ig	GenBank Accession Number:	BC018120	Purification Method:	Protein A purification
Size:	150ul	GeneID (NCBI):	5297	CloneNo.:	4E8C11
Source:	Mouse	UNIPROT ID:	P42356		
Isotype:	IgG2b	Full Name:	phosphatidylinositol 4-kinase, catalytic, alpha		
Immunogen Catalog Number:	AG17168	Calculated MW:	2044 aa, 231 kDa		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

PI4KA, also named as PIK4, PIK4CA, belongs to the PI3/PI4-kinase family and Type III PI4K subfamily. It acts on phosphatidylinositol (PtdIns) in the first committed step in the production of the second messenger inositol-1,4,5,-trisphosphate. PI4KA plays a role in the formation of membrane complexes where HCV replication takes place. (PMID:19605471).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

