

For Research Use Only

Albumin Monoclonal antibody

Catalog Number: 66051-2-Ig **1 Publications**



Basic Information

Catalog Number: 66051-2-Ig	GenBank Accession Number: BC034023	Purification Method: Protein G purification
Size: 150ul	GeneID (NCBI): 213	CloneNo.: 5E6B5
Source: Mouse	UNIPROT ID: P02768	
Isotype: IgG1	Full Name: albumin	
Immunogen Catalog Number: AG9885	Calculated MW: 609 aa, 69 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF

Species Specificity:
human

Cited Species:
human

Notable Publications

Author	Pubmed ID	Journal	Application
Zhizhao Deng	39680749	Adv Sci (Weinh)	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

