

For Research Use Only

COX7A2L Monoclonal antibody

Catalog Number: 66062-1-Ig 2 Publications



Basic Information

Catalog Number: 66062-1-Ig	GenBank Accession Number: BC007095	Purification Method: Protein G purification
Size: 150ul , Concentration: 760 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9167	CloneNo.: 3F9C2
Source: Mouse	UNIPROT ID: O14548	
Isotype: IgG1	Full Name: cytochrome c oxidase subunit VIIa polypeptide 2 like	
Immunogen Catalog Number: AG18768	Calculated MW: 114 aa, 13 kDa	
	Observed MW: 14 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Mina Davoudi	27997587	PLoS One	WB
Paula Garcia-Esparcia	27984680	Brain Pathol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

