

For Research Use Only

RBPJ Monoclonal antibody, PBS Only

Catalog Number: 66132-1-PBS

Featured Product



Basic Information

Catalog Number:

66132-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6302

GenBank Accession Number:

BC064976

GeneID (NCBI):

3516

UNIPROT ID:

Q06330

Full Name:

recombination signal binding protein for immunoglobulin kappa J region

Calculated MW:

56 kDa

Observed MW:

60 kDa

Purification Method:

Protein G purification

CloneNo.:

2F12G12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of RBPJ antibody (66132-1-Ig; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-RBPJ transfected HeLa cells. This data was developed using the same antibody clone with 66132-1-PBS in a different storage buffer formulation.



HepG2 cells were subjected to SDS PAGE followed by western blot with 66132-1-Ig (RBPJ antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66132-1-PBS in a different storage buffer formulation.