

For Research Use Only

IL-33 Monoclonal antibody, PBS Only

Catalog Number: 66235-2-PBS



Basic Information

Catalog Number:

66235-2-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

EG0313

GenBank Accession Number:

BC047085

GeneID (NCBI):

90865

UNIPROT ID:

O95760

Full Name:

interleukin 33

Calculated MW:

30 kDa

Purification Method:

Protein G purification

CloneNo.:

2B5B5

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

