

For Research Use Only

DIAPH3 Monoclonal antibody

Catalog Number: 66273-1-Ig



Basic Information

Catalog Number:

66273-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6631

GenBank Accession Number:

BC048963

GeneID (NCBI):

81624

UNIPROT ID:

Q9NSV4

Full Name:

diaphanous homolog 3 (Drosophila)

Calculated MW:

137 kDa

Purification Method:

Protein G purification

CloneNo.:

2F11G6

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

