

Phospho-TDP43 (Ser409/410) Monoclonal antibody, PBS Only

Catalog Number: 66318-1-PBS

Basic Information

Catalog Number: 66318-1-PBS	GenBank Accession Number: NM_007375	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 23435	CloneNo.: 1A2C1
Source: Mouse	UNIPROT ID: Q13148	
Isotype: IgG1	Full Name: TAR DNA binding protein	
	Calculated MW: 43 kDa	
	Observed MW: 45-50 kDa, 35 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse, rat

Background Information

The TARDBP gene encodes the TDP-43 protein, initially found to repress HIV-1 transcription by binding TAR DNA. TDP-43 has since been shown to bind RNA as well as DNA, and have multiple functions in transcriptional repression, translational regulation and pre-mRNA splicing. For instance, it is reported to regulate alternate splicing of the CTFR gene. In 2006 Neumann et al. found that hyperphosphorylated, ubiquitinated and/or cleaved forms of TDP-43, collectively known as pathological TDP-43, play a major role in the disease mechanisms of ubiquitin-positive, tau- and alpha-synuclein-negative frontotemporal dementia (FTLD-U) and in amyotrophic lateral sclerosis (ALS). Various forms of TDP-43 exist, including 18-35 kDa of cleaved C-terminal fragments, 45-50 kDa phospho-protein, 55 kDa glycosylated form, 75 kDa hyperphosphorylated form, and 90-300 kDa cross-linked form. (PMID: 17023659, 19823856, 21666678, 22193176). 66318-1-Ig is a phospho-dependent monoclonal antibody specifically recognizing phospho-TDP43 (409/410). It does not react with native TDP-43.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated and EA treated HeLa cells were subjected to SDS PAGE followed by western blot with 66318-1-Ig (Phospho-TDP43 (Ser409/410) Antibody) at dilution of 1:850 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66318-1-PBS in a different storage buffer formulation.