

For Research Use Only

Mouse IgG3 isotype control Monoclonal antibody



Catalog Number: 66360-4-Ig

Basic Information

Catalog Number: 66360-4-Ig	GenBank Accession Number: GeneID (NCBI):	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	Full Name:	CloneNo.: 1H4A5
Source: Mouse		Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG3		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB ,
Species Specificity: Megathura crenulata	

Background Information

This antibody is used as the isotype control of mouse IgG3. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

