

For Research Use Only

Mouse IgG3 isotype control Monoclonal antibody

Catalog Number: 66360-4-Ig



Basic Information

Catalog Number:

66360-4-Ig

GenBank Accession Number:

GeneID (NCBI):

Size:

150ul , Concentration: 1000 ug/ml by Nanodrop;

Full Name:

Source:

Mouse

Isotype:

IgG3

Purification Method:

Protein A purification

CloneNo.:

1H4A5

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

megathura crenulata

Positive Controls:

WB : ,

Background Information

This antibody is used as the isotype control of mouse IgG3. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66360-4-Ig (KLH antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.