

For Research Use Only

LCP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66465-2-PBS



Basic Information

Catalog Number:

66465-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3432

GenBank Accession Number:

BC016618

GeneID (NCBI):

3937

UNIPROT ID:

Q13094

Full Name:

Lymphocyte cytosolic protein 2 (SH2 domain containing leukocyte protein of 76kDa)

Calculated MW:

533 aa, 60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A12C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

66465-2-PBS targets LCP2 as part of a matched antibody pair:

MP51467-1: 66465-2-PBS capture and 66465-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

