

For Research Use Only

Rabbit Heavy chain of Rabbit IgG Monoclonal antibody, PBS Only



Catalog Number: 66467-1-PBS

Basic Information

Catalog Number:

66467-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

UNIPROT ID:

N/A

Full Name:

Observed MW:

50-55 kDa

Purification Method:

Protein G purification

CloneNo.:

1C9A5

Applications

Tested Applications:

WB, ELISA

Species Specificity:

rabbit

Background Information

This antibody specifically recognizes the heavy chain of Rabbit IgG.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rabbit serum (10nL) were subjected to SDS PAGE followed by western blot with 66467-1-Ig (mouse McAb to Heavy chain of Rabbit IgG) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66467-1-PBS in a different storage buffer formulation.