

For Research Use Only

HRPT2; CDC73 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66490-3-PBS



Basic Information

Catalog Number:

66490-3-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2959

GenBank Accession Number:

BC014351

GeneID (NCBI):

79577

UNIPROT ID:

Q6P1J9

Full Name:

cell division cycle 73, Paf1/RNA
polymerase II complex component,
homolog (S. cerevisiae)

Calculated MW:

531 aa, 61 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G4E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

66490-3-PBS targets HRPT2; CDC73 as part of a matched antibody pair:

MP51478-1: 66490-2-PBS capture and 66490-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51478-1, CDC73 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66490-2-PBS. Detection antibody: 66490-3-PBS.



Cytometric bead array standard curve of MP51478-1, CDC73 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66490-2-PBS. Detection antibody: 66490-3-PBS. Standard: Ag2959. Range: 3.125-400 ng/mL.