

For Research Use Only

SURVIVIN Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66495-4-PBS



Basic Information

Catalog Number: 66495-4-PBS	GenBank Accession Number: BC008718	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 332	CloneNo.: 3B9H7
Source: Mouse	UNIPROT ID: O15392	
Isotype: IgG1	Full Name: baculoviral IAP repeat-containing 5	
Immunogen Catalog Number: AG20958	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66495-4-PBS targets SURVIVIN as part of a matched antibody pair:

MP51208-2: 66495-4-PBS capture and 66495-3-PBS detection (validated in Cytometric bead array)

MP51208-4: 66495-4-PBS capture and 60760-1-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51208-2, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-4-PBS. Detection antibody: 66495-3-PBS. Standard:Ag20958. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51208-4, SURVIVIN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66495-4-PBS. Detection antibody: 60760-1-PBS. Standard:Ag20958. Range: 3.125-100 ng/mL