

For Research Use Only

SIAE Monoclonal antibody, PBS Only

Catalog Number: 66522-1-PBS



Basic Information

Catalog Number:

66522-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14624

GenBank Accession Number:

BC068450

GeneID (NCBI):

54414

UNIPROT ID:

Q9HAT2

Full Name:

sialic acid acetyltransferase

Calculated MW:

58 kDa

Observed MW:

70 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5F7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse , rat , pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66522-1-Ig (SIAE antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66522-1-PBS in a different storage buffer formulation.



pig colon tissue were subjected to SDS PAGE followed by western blot with 66522-1-Ig (SIAE antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66522-1-PBS in a different storage buffer formulation.