

For Research Use Only

# Ki-67 Monoclonal antibody

Catalog Number: 66555-1-Ig **3 Publications**



## Basic Information

<b>Catalog Number:</b> 66555-1-Ig	<b>GenBank Accession Number:</b> NM_002417	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 4288	<b>CloneNo.:</b> 4G8A6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P46013	
<b>Isotype:</b> IgG1	<b>Full Name:</b> antigen identified by monoclonal antibody Ki-67	
<b>Immunogen Catalog Number:</b> AG26266	<b>Calculated MW:</b> 359 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IHC, IF

**Species Specificity:**  
Human

**Cited Species:**  
human, rat

## Notable Publications

Author	Pubmed ID	Journal	Application
Xiaozhi Li	36316517	J Cancer Res Clin Oncol	IHC
Chaoqiong Ding	41747254	Cancer Discov	IF
Shaojun Fang	40368885	Light Sci Appl	IHC

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data