

For Research Use Only

MT-ND5 Monoclonal antibody, PBS Only

Catalog Number: 66613-1-PBS



Basic Information

Catalog Number:

66613-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

GenBank Accession Number:

AF465942

GeneID (NCBI):

4540

UNIPROT ID:

P03915

Full Name:

NADH dehydrogenase, subunit 5 (complex I)

Calculated MW:

67 kDa

Observed MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

1C12G8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66613-1-Ig (ND5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66613-1-PBS in a different storage buffer formulation.