

For Research Use Only

KIAA0284 Monoclonal antibody

Catalog Number: 66749-1-Ig



Basic Information

Catalog Number:

66749-1-Ig

GenBank Accession Number:

BC112928

Purification Method:

Protein A purification

Size:

150ul, Concentration: 1800 ug/ml by Nanodrop and 1000 ug/ml by Bradford method using BSA as the standard;

GeneID (NCBI):

283638

CloneNo.:

1D4B8

Source:

Mouse

UNIPROT ID:

Q9Y4F5

Recommended Dilutions:

WB 1:1000-1:4000

Isotype:

IgG2b

Full Name:

KIAA0284

Calculated MW:

1589 aa, 172 kDa

Immunogen Catalog Number:

AG26852

Observed MW:

172 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : HeLa cells, HEK-293 cells

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66749-1-Ig (KIAA0284 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.