

For Research Use Only

KIAA0284 Monoclonal antibody

Catalog Number: 66749-1-Ig



Basic Information

Catalog Number: 66749-1-Ig	GenBank Accession Number: BC112928	Purification Method: Protein A purification
Size: 150ul , Concentration: 1800 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 283638	CloneNo.: 1D4B8
Source: Mouse	Full Name: KIAA0284	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG2b	Calculated MW: 1589 aa, 172 kDa	
Immunogen Catalog Number: AG26852	Observed MW: 172 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66749-1-Ig (KIAA0284 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.