

For Research Use Only

S100A8 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66853-2-PBS



Basic Information

Catalog Number: 66853-2-PBS	GenBank Accession Number: BC005928	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6279	CloneNo.: 1B11F6
Source: Mouse	UNIPROT ID: P05109	
Isotype: IgG1	Full Name: S100 calcium binding protein A8	
Immunogen Catalog Number: AG8517	Calculated MW: 93 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

66853-2-PBS targets S100A8 as part of a matched antibody pair:

MP80012-1: 82531-3-PBS capture and 66853-2-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

