

For Research Use Only

TRIM45 Monoclonal antibody, PBS Only



Catalog Number: 66896-1-PBS

Basic Information

Catalog Number:

66896-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28300

GenBank Accession Number:

BC034943

GeneID (NCBI):

80263

UNIPROT ID:

Q9H8W5

Full Name:

tripartite motif-containing 45

Calculated MW:

562aa, 62 kDa; 580aa, 64 kDa

Observed MW:

62-64 kDa

Purification Method:

Protein G purification

CloneNo.:

3A7H11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 66896-1-Ig (TRIM45 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66896-1-PBS in a different storage buffer formulation.