

For Research Use Only

IKZF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66966-2-PBS



Basic Information

Catalog Number: 66966-2-PBS	GenBank Accession Number: BC018349	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10320	CloneNo.: 1C1G12
Source: Mouse	UNIPROT ID: Q13422	
Isotype: IgG1	Full Name: IKAROS family zinc finger 1 (Ikaros)	
Immunogen Catalog Number: AG28637	Calculated MW: 477 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

66966-2-PBS targets IKZF1 as part of a matched antibody pair:

MP51595-1: 66966-2-PBS capture and 66966-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS.



Cytometric bead array standard curve of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS. Standard: Ag28637. Range: 1.563-200 ng/mL.