

For Research Use Only

EDN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67008-4-PBS



Basic Information

Catalog Number: 67008-4-PBS	GenBank Accession Number: BC009720	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1000 µg/ml by 1906 Nanodrop;	GeneID (NCBI): UNIPROT ID: P05305	CloneNo.: 3D10B8
Source: Mouse	Full Name: endothelin 1	
Isotype: IgG2a	Calculated MW: 212 aa, 24 kDa	
Immunogen Catalog Number: AG28334		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

67008-4-PBS targets EDN1 as part of a matched antibody pair:

MP51065-2: 67008-4-PBS capture and 67008-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51065-2, EDN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67008-4-PBS. Detection antibody: 67008-5-PBS.



Cytometric bead array standard curve of MP51065-2, EDN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67008-4-PBS. Detection antibody: 67008-5-PBS. Standard:Ag28334. Range: 1.563-200 ng/mL.