

For Research Use Only

LSM12 Monoclonal antibody, PBS Only



Catalog Number: 67029-1-PBS

Featured Product

Basic Information

Catalog Number:

67029-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28782

GenBank Accession Number:

BC044587

GeneID (NCBI):

124801

UNIPROT ID:

Q3MHD2

Full Name:

LSM12 homolog (S. cerevisiae)

Observed MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

2E3G3

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

Human, Mouse, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67029-1-Ig (LSM12 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67029-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67029-1-Ig (LSM12 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67029-1-PBS in a different storage buffer formulation.