

For Research Use Only

# NGF Monoclonal antibody

Catalog Number: 67036-1-Ig



## Basic Information

Catalog Number:

67036-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14564

GenBank Accession Number:

BC032517

GeneID (NCBI):

4803

UNIPROT ID:

P01138

Full Name:

nerve growth factor (beta polypeptide)

Calculated MW:

241 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

1C3B11

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

