

For Research Use Only

DDR2 Monoclonal antibody, PBS Only

Catalog Number: 67126-1-PBS



Basic Information

Catalog Number: 67126-1-PBS	GenBank Accession Number: NM_001014796	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4921	CloneNo.: 1E9A10
Source: Mouse	UNIPROT ID: Q16832	
Isotype: IgG1	Full Name: discoidin domain receptor tyrosine kinase 2	
Immunogen Catalog Number: AG28340	Calculated MW: 97 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 67126-1-Ig (DDR2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67126-1-PBS in a different storage buffer formulation.