

For Research Use Only

# PLVAP Monoclonal antibody, PBS Only



Catalog Number: 67158-1-PBS

## Basic Information

Catalog Number:

67158-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6359

GenBank Accession Number:

BC056414

GeneID (NCBI):

83483

UNIPROT ID:

Q9BX97

Full Name:

plasmalemma vesicle associated protein

Calculated MW:

51 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

1D4A6

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

Human, Pig

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 67158-1-Ig (PLVAP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67158-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67158-1-Ig (PLVAP antibody) at dilution of 1:4000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67158-1-PBS in a different storage buffer formulation.