

For Research Use Only

# DRG2 Monoclonal antibody, PBS Only



Catalog Number: 67164-1-PBS

Featured Product

## Basic Information

Catalog Number:

67164-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6878

GenBank Accession Number:

BC000493

GeneID (NCBI):

1819

UNIPROT ID:

P55039

Full Name:

developmentally regulated GTP binding protein 2

Calculated MW:

41 kDa

Observed MW:

41 kDa

Purification Method:

Protein A purification

CloneNo.:

3A3E1

## Applications

Tested Applications:

WB, IF, Indirect ELISA

Species Specificity:

Human, mouse, rat, pig

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67164-1-Ig (DRG2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67164-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed U2OS cells using DRG2 antibody (67164-1-Ig, Clone: 3A3E1) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 67164-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed Saos-2 cells using DRG2 antibody (67164-1-Ig, Clone: 3A3E1) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67164-1-PBS in a different storage buffer formulation.