

For Research Use Only

KIF20A Monoclonal antibody, PBS Only

Catalog Number: 67190-1-PBS



Basic Information

Catalog Number: 67190-1-PBS	GenBank Accession Number: BC012999	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 10112	CloneNo.: 2A8E3
Source: Mouse	UNIPROT ID: O95235	
Isotype: IgG1	Full Name: kinesin family member 20A	
Immunogen Catalog Number: AG8916	Calculated MW: 890 aa, 100 kDa	
	Observed MW: 100 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
Human, Rat, Mouse

Background Information

KIF20A, also known as RAB6KIFL or MKlp2, is a kinesin originally identified as a binding partner for the rab6 GTPase involved in protein transport at the Golgi apparatus. KIF20A belongs to a large family of motor proteins that accumulate in mitotic cells and is required for normal cell division. KIF20A is thought to be involved in carcinogenesis. In normal tissues KIF20A is almost undetectable while it is overexpressed in various cancers. KIF20A has been identified as a novel promising candidate for anticancer immunotherapy.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only




For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67190-1-Ig (KIF20A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 67190-1-PBS in a different storage buffer formulation.



WB result of KIF20A antibody (67190-1-Ig; 1:10000; incubated at room temperature for 1.5 hours) with sh-Control and sh-KIF20A transfected HeLa cells. This data was developed using the same antibody clone with 67190-1-PBS in a different storage buffer formulation.