

For Research Use Only

DTX2 Monoclonal antibody

Catalog Number: 67209-1-Ig



Basic Information

Catalog Number:

67209-1-Ig

GenBank Accession Number:

BC064535

Purification Method:

Protein A purification

Size:

150ul, Concentration: 2400 ug/ml by Nanodrop and 1000 ug/ml by Bradford method using BSA as the standard;

GeneID (NCBI):

113878

CloneNo.:

2D4F8

UNIPROT ID:

Q86UW9

Recommended Dilutions:

WB 1:1000-1:6000

Source:

Mouse

Full Name:

deltex homolog 2 (Drosophila)

Isotype:

IgG2a

Calculated MW:

253aa, 28 kDa; 575aa, 64 kDa

Immunogen Catalog Number:

AG12022

Observed MW:

68 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells, A549 cells

Species Specificity:

Human

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67209-1-Ig (DTX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.