

For Research Use Only

# MMP14 Monoclonal antibody

Catalog Number: 67225-1-Ig



## Basic Information

Catalog Number:

67225-1-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6062

GenBank Accession Number:

BC064803

GeneID (NCBI):

4323

UNIPROT ID:

P50281

Full Name:

matrix metalloproteinase 14  
(membrane-inserted)

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

2F11A11

## Applications

Tested Applications:

ELISA

Species Specificity:

Human, mouse, rat

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

